

Your patient's smile is as special as they are. Why hide it behind wires? Thanks to Clearbow™ labial bows, you no longer need to. NorthStar is one of the nation's leading suppliers of Clearbow™ retainers. They are the perfect solution for an aesthetically pleasing, natural looking smile, while utilizing the benefits of a traditional retainer.

Aesthetic labial bows are not a new invention. Perhaps you've tried them in the past, only to experience the following: broken bows, patients complaining of stains, difficulty in adapting the bow to your patient's teeth, and high costs. Thanks to Clearbow™ those days are over.

The Basics

Clearbow™ is fabricated in Canada using food-grade materials made only in the United States. The bow material is made of polyethylene terephthalate (PET), which is the same BPA-free plastic used to make soda and juice bottles. The wire used for the Coiltight-Joints™ is grade 316 stainless steel.



The Benefits

Clearbow™ provides you and your patients with many benefits, including:

- Unparalleled aesthetics for the wearing pleasure of your patient.
- 3-prong plier adjustable, so you can adapt it the same as you would a steel labial bow.
- 2.75mm of full labial contact and a concave bow structure provide the patient with optimum tooth contact for superior retention.
- Unaffected by boiling water, and resistant to temperatures up to 250° F.
- Gentle on teeth. Will not stain or scratch tooth enamels, restorations, or veneers.
- Virtually resistant to stains. Patients can enjoy beverages like coffee, red wine, and cranberry juice without worry.
- Affordable. Clearbow™ is available for a fraction of the price of other aesthetic labial bows.

The Guarantee

Clearbow's Coiltight-Joints™ can withstand a combined static weight of 55 lbs, and the plastic strap is virtually indestructible. That's why the makers of Clearbow™ guarantee that the plastic strap of the ClearBow™ assembly not to break or detach from the Coiltight-Joints™ when the orthodontic retainer is used according to the doctor's instructions.



To learn more about Clearbow™ please visit their website @ <http://www.clearbow.com>

Sources: Clearbow. (n.d.). Retrieved April 9, 2015, from <http://www.clearbow.com/doctors/>
Advantages. (n.d.). Retrieved April 9, 2015, from <http://www.clearbow.com/doctors/>

Clearbow is a trademark of PNG Orthodontic Specialties, Ltd.
Coiltight-Joints is a trademark of PNG Orthodontic Specialties, Ltd.