

# Bite Registrations and Models Required for Full Warranty

<u>APPLIANCE</u>	<u>ONE ARCH</u>	<u>BOTH ARCHES W/WAX BITE</u>	<u>A/P RELATION</u>	<u>OPENING</u>
<b>Retention</b>				
Standard Designs	Yes	No	N/A	N/A
<b>Finishing and Aligning</b>				
Spring Aligner (Standard Designs)		Yes	Centric	None
Star Aligner		Yes	Centric	None
Jump Start Aligner		Yes	Centric	None
<b>Health and Aesthetics</b>				
Ozone Therapy Tray	Yes	No	N/A	N/A
<b>Safety</b>				
Mouthguard (Single Arch)	Yes	No	N/A	N/A
<b>Arch Development</b>				
Rapid Palatal Expander	Yes	No	N/A	N/A
Rapid Palatal Expander (Bondable)		No*	Centric	2-3mm Posterior
3D™ Designs	Yes	No	N/A	N/A
Haas Expander	Yes	No	N/A	N/A
Transverse (Standard Designs)		Yes	Centric	2-3mm Posterior
Sagittal (Standard Designs)		Yes	Centric	3mm Posterior
Schwarz	Yes	No	N/A	N/A
<b>Functional Therapy</b>				
Activator (Standard Designs)		Yes	Edge-To-Edge	6-9mm Posterior
Bionator I		Yes	Edge-To-Edge	2-3mm Anterior
Bionator II & to Maintain		Yes	Edge-To-Edge	2-3mm Posterior
Herbst®		Yes	As Needed	No Additional
Frankel I & II		Yes	3mm	3mm Posterior
Frankel III		Yes	Edge-To-Edge	2-3mm Posterior
Twin Block		Yes	Edge-To-Edge	3-4mm Posterior
<b>Habit Control</b>				
Standard Designs		Yes	Centric	None
<b>Splint Therapy</b>				
Standard Designs		Yes	As Needed	Less than 3mm Posterior
<b>Design Additions</b>				
Pontic		No*	Centric	None
Edentulous Block		No*	Centric	None
Posterior Coverage		No*	Centric	None
Tongue Guard		Yes	Centric	None
Resets		Yes	Centric	None
Flange		Yes	Centric	None
Headgear Tubes		Yes	Centric	2-3mm Posterior

\*Although not necessary for fabrication, opposing models with a wax bite may produce a better fitting appliance.

