

DIGITAL SERVICES

DIGITAL STUDY MODELS

FROM ALGINATE IMPRESSIONS

STANDARD	40
PROGRESS	35
UNBASED	29

FROM HARD MODELS

STANDARD	26
PROGRESS	21
UNBASED	15

FROM DIGITAL IMPRESSIONS

STANDARD	16
PROGRESS	16
UNBASED	16

DSM REPORTS

STANDARD (choose .jpg or .pdf)	4
OB/OJ (choose .jpg or .pdf)	5

3D MODEL PRINTING

FOR NS APPLIANCE CONSTRUCTION	FREE
RESIN MODEL PRINTING (per gram)	CALL

FINISHING & ALIGNING

ANTERIOR SPRING	60
STANDARD SPRING COMPLETE	86
INVISIBLE ALIGNER W/RESETS	72

VIRTUAL DIAGNOSTIC SET-UP

SINGLE ARCH	63*
DUAL ARCH	80*

*VIRTUAL DIAGNOSTIC FEE IS WAIVED WITH PURCHASE OF A SET OF ALIGNERS

SERVICES & FEES

BAND/BRACKET REMOVAL (per arch)	5
SINGLE ARCH SCAN	6
PLASTER STUDY MODEL FROM DSM	113
PLACE DR. PROVIDED BANDS (ea.)	8
PLACE NORTHSTAR BANDS (ea.)	12
CONVERT FILE TO .STL (ea.)	7
DIGITAL CANCELLATION FEE (per case)	6.5
ORTHOVIEWER™ SOFTWARE	FREE
BOXES FOR BAND SUBMISSION	FREE

ADDITIONAL INFORMATION

Supported Scanners

NorthStar currently accepts scans in .stl format from the following intraoral scanners: iTero® family of scanners, 3Shape Trios™, Carestream CS 3500™, 3M True Definition™, Cerec Ortho™, Lythos™ and all other popular brands.

Scan Recommendations

In order to ensure the best results, your submitted scans should capture as much detail from the arch and surrounding tissues as is reasonably possible. NorthStar requires all submitted scans to meet the following minimum quality requirements:

Upper arch - Scan full upper arch and palate extending mesial the 6s. Must extend at least 5mm below the cervical line.

Lower arch - Scan entire arch to floor of mouth and back to most distal molar. Must extend at least 5mm below the cervical line.

Please contact our Digital Services department at 800.346.0011 for additional assistance.



218 Industrial Park Rd.
Park Rapids, MN 56470
800.346.0011
www.northstarorthodontics.com
info@northstarorthodontics.com